BioPak 240 R BENCH CONTESTANT					
	VISUAL APPARATUS CHECKS		TEST APPARATUS		
1	/	neck if ok	1	C	heck if ok
U	UF	PPER HOUSING			
/	LC	WER HOUSING	V	Ze	ero Adjust the Mag. Gauge
V	CE	NTER SECTION ASSEMBLY		/	DNSTANT FLOW TEST
/	CA	RBON DIOXIDE SCRUBBER		1	State Reading per Rev.D p. 10
/	SIN	JULATE ADDING SCRUBBERS IF NOT PRESENT		/ DE	EMAND VALVE TEST Flow stops 1-3 Sec.
/	CE	NTER SECTION LID ASSEMBLY		1	AFROTAIN AND AND AND AND AND AND AND AND AND AN
<b>V</b>	SIN	NULATE ADDING ICE IF NOT PRESENT	U	1	INT VALVE TEST Flow no more 1-2 Sec.
1	СУ	LINDER TEST		1	State Reading 2 inches of water or less
	но	SES	1	LO	W PRESSURE LEAK TEST No more than02" WC
J	FAC	CEPIECE TEST	J	ſ	S GAUGE AND TRIM TEST
			V		State Reading Observe lights/gauges +/- 10%
			_	Ніс	SH PRESSURE LEAK TEST
L,	COI	NNECTIONS CONNECTED OR GOOD	V	,	W PRESSURE ALARM TEST
1	Ven	t Valve Assembly - Connected or good	V	/	State Reading Alarm 650-1000 psig
0	Diap	ohragm Worm Gear - Connected or good			Alaim 050-1000 psig
V	Flov	v Restrictor - Connected or good			
<b>S</b>	Brea	athing Hose Worm Gear - Connected or good			
<b>V</b>	Add / Constant Fittings - Connected or good				
5	Center Section Lid - Connected or good				
1	Center Section Push Pins - Connected or good				
J	Cylinder Connection - Connected or good				
	Adapter to Facepiece - Connected or good				
J	Test	Fixture Connections - Connected or good			
VI	VC	Work			
	V	VISUAL	VI	VC	TEST
-	1	Missing user right OTS Mounting Plate Missing PCM			<del>Demand Valve te</del> ak
	V			V	Gage Disconnected
	V /	Missing Flow Restrictor Assembly (Inc Washer)  Missing Vent Spacer		V.	Missing Center Section O-ring/Low Pressure leak
	V			$\checkmark$	Missing Hose Adapter O-ring/Low Pressure leak
	V	Missing Hose Adapter O-ring			
,		Bottle Hydro test is out of date/Must remove bodok seal for new bottle			